

User Manual

Doc No: 910-0040 Rev A Dt.01/20/02

Model No: *GTM3TXX-YY-ZZZZ* (*XX=35/41/48/57*, *YY = 3 to 48*, *ZZZZ = PER CERTIFICATE*)

This product is a component for use with Medical Equipment. It is to be used with equipment for which it was designed to be used with.

WARNING: *Do not use this power supply in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.*

A. Introduction:

Model *GTM3TXX-YY-ZZZZ* is a Linear Power Supply, intended to operate from 220/230/240 VAC AT 50 Hz . In model number *XX* indicates case size *YY* indicate output voltage and *ZZZZ* indicate output current in milli-Amperes.

B. Setup and Operation:

1. Connect the power supply's output connector to power inlet of intended equipment.
2. Plug the power supply (Or Plug of Input Cord) into appropriate power socket.
3. power Supply will be immediately operational.

C. Specifications:

Please refer to product specification for product details.

CAUTION:

1. *Power Supply should not be operated outside the specifications.*
2. *Do Not use the power supply in wet locations.*
3. *Do Not open the power supply. Hazardous voltage exists inside the power supply.*